

RANBP17 Antibody

Catalog No: #34957

Package Size: #34957-1 50ul #34957-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RANBP17 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total RANBP17 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from N-terminal of human RANBP17.
Target Name	RANBP17
Other Names	Ran-binding protein 17; RANBP17;
Accession No.	Swiss-Prot: Q9H2T7NCBI Gene ID: 64901
Uniprot	Q9H2T7
GeneID	64901;
SDS-PAGE MW	125kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from Jurkat cells, using RANBP17 antibody #34957.

Background

May function as a nuclear transport receptor By similarity.

Koch P., Biochem. Biophys. Res. Commun. 278:241-249(2000).

Kutay U., J. Biol. Chem. 275:40163-40168(2000).

Ota T., Nat. Genet. 36:40-45(2004).

Note: This product is for in vitro research use only