

RANBP6 Antibody

Catalog No: #34958

Package Size: #34958-1 50ul #34958-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RANBP6 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total RANBP6 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human RANBP6.
Target Name	RANBP6
Other Names	RAN binding 6; RBNP6;
Accession No.	Swiss-Prot: O60518NCBI Gene ID: 26953
Uniprot	O60518
GeneID	26953;
SDS-PAGE MW	125kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from K562 cells, COLO cells and COS-7 cells, using RANBP6 antibody #34958.

Background

May function in nuclear protein import as nuclear transport receptor.

Goerlich D., Submitted (DEC-1997) to the EMBL/GenBank/DDBJ databases.

The German cDNA consortium; Submitted (JUN-2003) to the EMBL/GenBank/DDBJ databases.

Humphray S.J., Nature 429:369-374(2004).

Note: This product is for in vitro research use only