

G3BP2 Antibody

Catalog No: #34963

Package Size: #34963-1 50ul #34963-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	G3BP2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total G3BP2 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from C-terminal of human G3BP2.
Target Name	G3BP2
Other Names	G3B2; G3BP2; GAP SH3-domain binding protein 2; KIAA0660; Ras-GTPase-activating protein binding protein 2
Accession No.	Swiss-Prot: Q9UN86NCBI Gene ID: 9908
Uniprot	Q9UN86
GeneID	9908;
SDS-PAGE MW	55kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from A549 cells, HUVEC cells and HeLa cells, using G3BP2 antibody #34963.

Background

Probable scaffold protein that may be involved in mRNA transport Potential.

Kennedy D., Submitted (APR-1999) to the EMBL/GenBank/DDBJ databases.

Ishikawa K., DNA Res. 5:169-176(1998).

Kennedy D., Submitted (MAR-1998) to the EMBL/GenBank/DDBJ databases.

Note: This product is for in vitro research use only