

RAB41 Antibody

Catalog No: #34971

Package Size: #34971-1 50ul #34971-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RAB41 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total RAB41 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human RAB41.
Target Name	RAB41
Other Names	RAB41; member RAS homolog family;
Accession No.	Swiss-Prot: Q5JT25NCBI Gene ID: 347517
Uniprot	Q5JT25
GeneID	347517;
SDS-PAGE MW	24kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:50~1:100

Images



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue using RAB41 antibody #34971.

Background

Required for normal Golgi ribbon organization and ER-to-Golgi trafficking.

Ross M.T., Nature 434:325-337(2005).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

Note: This product is for in vitro research use only