

RAB6A Antibody

Catalog No: #34972

Package Size: #34972-1 50ul #34972-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RAB6A Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total RAB6A protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human.
Target Name	RAB6A
Other Names	gtp-binding protein rab6c; rab gtpase rab6a; rab6; rab6a; member ras oncogene family
Accession No.	Swiss-Prot: P20340NCBI Gene ID: 5870
Uniprot	P20340
GeneID	5870;
SDS-PAGE MW	24kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from A549 cells, using RAB6A antibody #34972.

Background

Protein transport. Regulator of membrane traffic from the Golgi apparatus towards the endoplasmic reticulum (ER). Has a low GTPase activity.

Zahraoui A., J. Biol. Chem. 264:12394-12401(1989).

Caillol N., J. Cell. Biochem. 79:628-647(2000).

Shan J., Gene 257:67-75(2000).

Note: This product is for in vitro research use only