

RAB7L1 Antibody

Catalog No: #34974

Package Size: #34974-1 50ul #34974-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RAB7L1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total RAB7L1 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human RAB7L1.
Target Name	RAB7L1
Other Names	RAB7; member RAS oncogene family-like 1; RAB7L;
Accession No.	Swiss-Prot: O14966NCBI Gene ID: 8934
Uniprot	O14966
GeneID	8934;
SDS-PAGE MW	28kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from HT-29 cells, using RAB7L1 antibody #34974.

Background

Shimizu F., Cytogenet. Cell Genet. 77:261-263(1997).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

Note: This product is for in vitro research use only