

RREB1 Antibody

Catalog No: #34975

Package Size: #34975-1 50ul #34975-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RREB1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total RREB1 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human RREB1.
Target Name	RREB1
Other Names	RREB-1; Raf-responsive zinc finger protein LZ321; Zinc finger motif-enhancer binding-protein 1; Zep-1; Finger protein in nuclear bodies
Accession No.	Swiss-Prot: Q92766NCBI Gene ID: 6239
Uniprot	Q92766
GeneID	6239;
SDS-PAGE MW	180kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:50~1:100

Images



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue using RREB1 antibody #34975.

Background

Transcription factor that binds specifically to the RAS-responsive elements (RRE) of gene promoters. May be involved in Ras/Raf-mediated cell differentiation by enhancing calcitonin expression. Represses the angiotensinogen gene. Negatively regulates the transcriptional activity of AR. Potentiates the transcriptional activity of NEUROD1.

Fujimoto-Nishiyama A., Gene 195:267-275(1997).

Date S., Int. J. Mol. Med. 13:637-642(2004).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

Note: This product is for in vitro research use only