

RFPL4A Antibody

Catalog No: #34985

Package Size: #34985-1 50ul #34985-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RFPL4A Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total RFPL4A protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human RFPL4A.
Target Name	RFPL4A
Other Names	Ret finger protein-like 4A; RING finger protein 210; RFPL4A; RFPL4; RNF210
Accession No.	Swiss-Prot: A6NLU0NCBI Gene ID: 342931
Uniprot	A6NLU0
GeneID	342931;
SDS-PAGE MW	42kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from HUVEC cells, using RFPL4A antibody #34985.

Background

Grimwood J., Nature 428:529-535(2004).

Gruetter C., FEBS Lett. 580:99-106(2006).

Note: This product is for in vitro research use only