

RHG11A Antibody

Catalog No: #34987

Package Size: #34987-1 50ul #34987-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RHG11A Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total RHG11A protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human RHG11A.
Target Name	RHG11A
Other Names	Rho-type GTPase-activating protein 11A; GAP (1-12); KIAA0013; MGC70740; RHGBA
Accession No.	Swiss-Prot: Q6P4F7NCBI Gene ID: 9824
Uniprot	Q6P4F7
GeneID	9824;
SDS-PAGE MW	110kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from NIH/3T3 cells, using RHG11A antibody #34987.

Background

Nomura N., DNA Res. 1:27-35(1994).

Nousiainen M., Proc. Natl. Acad. Sci. U.S.A. 103:5391-5396(2006).

Wang B., Anal. Chem. 80:9526-9533(2008)

Note: This product is for in vitro research use only