

ARHGEF19 Antibody

Catalog No: #34994

Package Size: #34994-1 50ul #34994-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ARHGEF19 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total ARHGEF19 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human ARHGEF19.
Target Name	ARHGEF19
Other Names	Rho guanine nucleotide exchange factor 19; ARHGEF19;
Accession No.	Swiss-Prot: Q8IW93NCBI Gene ID: 128272
Uniprot	Q8IW93
GeneID	128272;
SDS-PAGE MW	90kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:50~1:100

Images



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue using ARHGEF19 antibody #34994.

Background

Acts as guanine nucleotide exchange factor (GEF) for RhoA GTPase By similarity.

Ota T., Nat. Genet. 36:40-45(2004).

Gregory S.G., Nature 441:315-321(2006).

Note: This product is for in vitro research use only