

RIMKA Antibody

Catalog No: #35001

Package Size: #35001-1 50ul #35001-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RIMKA Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total RIMKA protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human RIMKA.
Target Name	RIMKA
Other Names	ribosomal protein S6 modification-like A; RIMKA;
Accession No.	Swiss-Prot: Q8IXN7NCBI Gene ID: 284716
Uniprot	Q8IXN7
GeneID	284716;
SDS-PAGE MW	43kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from Jurkat cells, HuvEc cells and COS cells, using RIMKA antibody #35001.

Background

Catalyzes the synthesis of N-acetyl-L-aspartyl-L-glutamate (NAAG) and N-acetyl-L-aspartyl-L-glutamyl-L-glutamate.

Gregory S.G., Nature 441:315-321(2006).

The MGC Project Team; Genome Res. 14:2121-2127(2004)

Note: This product is for in vitro research use only