

RRS1 Antibody

Catalog No: #35002

Package Size: #35002-1 50ul #35002-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RRS1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous levels of total RRS1 protein
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human RRS1.
Target Name	RRS1
Other Names	ribosome biogenesis regulatory protein; RRR;
Accession No.	Swiss-Prot: Q15050NCBI Gene ID: 23212
Uniprot	Q15050
GeneID	23212;
SDS-PAGE MW	41kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from 293 cells and HuvEc cells, using RRS1 antibody #35002.

Background

Involved in ribosome biogenesis By similarity.

Nagase T., DNA Res. 2:37-43(1995).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

Bienvenut W.V.; Submitted (AUG-2005) to UniProtKB.

Note: This product is for in vitro research use only