

SCARF2 Antibody

Catalog No: #35010

Package Size: #35010-1 50ul #35010-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SCARF2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total SCARF2 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human SCARF2.
Target Name	SCARF2
Other Names	Scavenger receptor class F member 2 Scavenger receptor expressed by endothelial cells 2 protein; SREC-II; SRECRP-1 SCARF2;
Accession No.	Swiss-Prot: Q96GP6NCBI Gene ID: 91179
Uniprot	Q96GP6
GeneID	91179;
SDS-PAGE MW	92kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:50~1:100

Images



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue using SCARF2 antibody #35010.

Background

Probable adhesion protein, which mediates homophilic and heterophilic interactions. In contrast to SCARF1, it poorly mediates the binding and degradation of acetylated low density lipoprotein (Ac-LDL) By similarity.

Ishii J., J. Biol. Chem. 277:39696-39702(2002).

Sjoeblom T., Science 314:268-274(2006).

Note: This product is for in vitro research use only