

SERGEF Antibody

Catalog No: #35012

Package Size: #35012-1 50ul #35012-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SERGEF Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total SERGEF protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human SERGEF.
Target Name	SERGEF
Other Names	Guanine nucleotide exchange factor-related protein; Deafness locus-associated putative guanine nucleotide exchange factor; DelGEF;
Accession No.	Swiss-Prot: Q9UGK8NCBI Gene ID: 26297
Uniprot	Q9UGK8
GeneID	26297;
SDS-PAGE MW	50kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from NIH/3T3 cells, using SERGEF antibody #35012.

Background

Probable guanine nucleotide exchange factor (GEF), which may be involved in the secretion process.

Uhlmann J., FEBS Lett. 460:153-160(1999).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

Sjoelinder M., FEBS Lett. 532:211-215(2002).

Note: This product is for in vitro research use only