

SPINK6 Antibody

Catalog No: #35018

Package Size: #35018-1 50ul #35018-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SPINK6 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total SPINK6 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human SPINK6.
Target Name	SPINK6
Other Names	Serine protease inhibitor Kazal-type 6; SPINK6;
Accession No.	Swiss-Prot: Q6UWN8NCBI Gene ID: 404203
Uniprot	Q6UWN8
GeneID	404203;
SDS-PAGE MW	9kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:50~1:100

Images



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue using SPINK6 antibody #35018.

Background

Serine protease inhibitor selective for kallikreins. Efficiently inhibits KLK4, KLK5, KLK6, KLK7, KLK12, KLK13 and KLK14. Doesn't inhibit KLK8.

Clark H.F., Genome Res. 13:2265-2270(2003).

Zhang Z.,Protein Sci. 13:2819-2824(2004).

Note: This product is for in vitro research use only