

PPP2R3C Antibody

Catalog No: #35033

Package Size: #35033-1 50ul #35033-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	PPP2R3C Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total PPP2R3C protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human PPP2R3C.
Target Name	PPP2R3C
Other Names	Protein phosphatase subunit G5PR; Rhabdomyosarcoma antigen MU-RMS-40.6A; Rhabdomyosarcoma antigen Mu-RMS-40.6C;
Accession No.	Swiss-Prot: Q969Q6NCBI Gene ID: 55012
Uniprot	Q969Q6
GeneID	55012;
SDS-PAGE MW	53kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from HuvEc cells, using PPP2R3C antibody #35033.

Background

May regulate MCM3AP phosphorylation through phosphatase recruitment. May play a role in the activation-induced cell death of B-cells By similarity.

Kamnasaran D., Genomics 85:608-621(2005).

Ota T., Nat. Genet. 36:40-45(2004).

Mural R.J., Submitted (SEP-2005) to the EMBL/GenBank/DDBJ databases.

Note: This product is for in vitro research use only