

SH3TC1 Antibody

Catalog No: #35044

Package Size: #35044-1 50ul #35044-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SH3TC1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total SH3TC1 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human SH3TC1.
Target Name	SH3TC1
Other Names	SH3 domain and tetratricopeptide repeats-containing protein 1; SH3TC1;
Accession No.	Swiss-Prot: Q8TE82NCBI Gene ID: 54436
Uniprot	Q8TE82
GeneID	54436;
SDS-PAGE MW	147kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from K562 cells, using SH3TC1 antibody #35044.

Background

Ota T., Nat. Genet. 36:40-45(2004).

Hillier L.W., Nature 434:724-731(2005).

Mayya V., Sci. Signal. 2:RA46-RA46(2009).

Note: This product is for in vitro research use only