

SH3GL2 Antibody

Catalog No: #35047

Package Size: #35047-1 50ul #35047-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SH3GL2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total SH3GL2 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human SH3GL2.
Target Name	SH3GL2
Other Names	Endophilin-1; SH3 domain-containing GRB2-like protein 2; SH3 domain protein 2A; EEN-B1;
Accession No.	Swiss-Prot: Q99962NCBI Gene ID: 6456
Uniprot	Q99962
GeneID	6456;
SDS-PAGE MW	39kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from K562 cells and HeLa cells, using SH3GL2 antibody #35047.

Background

Implicated in synaptic vesicle endocytosis. May recruit other proteins to membranes with high curvature.

Giachino C., Genomics 41:427-434(1997).

So C.W., Submitted (MAR-1998) to the EMBL/GenBank/DDBJ databases.

Halleck A., Submitted (JUN-2004) to the EMBL/GenBank/DDBJ databases.

Note: This product is for in vitro research use only