

SLC5A6 Antibody

Catalog No: #35064

Package Size: #35064-1 50ul #35064-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SLC5A6 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total SLC5A6 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human SLC5A6.
Target Name	SLC5A6
Other Names	SMVT; solute carrier family 5 (sodium-dependent vitamin transporter); member 6;
Accession No.	Swiss-Prot: Q9Y289NCBI Gene ID: 8884
Uniprot	Q9Y289
GeneID	8884;
SDS-PAGE MW	68kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from HepG2 cells and HUVEC cells, using SLC5A6 antibody #35064.

Background

Transports pantothenate, biotin and lipoate in the presence of sodium.

Wang H., J. Biol. Chem. 274:14875-14883(1999).

Prasad P.D., Arch. Biochem. Biophys. 366:95-106(1999).

Note: This product is for in vitro research use only