

SPTBN5 Antibody

Catalog No: #35072

Package Size: #35072-1 50ul #35072-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SPTBN5 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total SPTBN5 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human SPTBN5.
Target Name	SPTBN5
Other Names	Spectrin beta chain; brain 4; Spectrin; non-erythroid beta chain 4; Beta-V spectrin
Accession No.	Swiss-Prot: Q9NRC6NCBI Gene ID: 51332
Uniprot	Q9NRC6
GeneID	51332;
SDS-PAGE MW	417kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:50~1:100

Images



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue using SPTBN5 antibody #35072.

Background

Belongs to the spectrin family.

Stabach P.R., J. Biol. Chem. 275:21385-21395(2000).

Note: This product is for in vitro research use only