

## SFRS11 Antibody

Catalog No: #35073

Package Size: #35073-1 50ul #35073-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

|                       |                                                                                                                                                    |
|-----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | SFRS11 Antibody                                                                                                                                    |
| Host Species          | Rabbit                                                                                                                                             |
| Clonality             | Polyclonal                                                                                                                                         |
| Purification          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                              |
| Applications          | WB                                                                                                                                                 |
| Species Reactivity    | Hu Ms                                                                                                                                              |
| Specificity           | The antibody detects endogenous levels of total SFRS11 protein.                                                                                    |
| Immunogen Type        | Peptide                                                                                                                                            |
| Immunogen Description | Synthesized peptide derived from internal of human SFRS11.                                                                                         |
| Target Name           | SFRS11                                                                                                                                             |
| Other Names           | Arginine-rich 54 kDa nuclear protein; p54; SFR11; SFRS11; Splicing factor arginine/serine-rich 11                                                  |
| Accession No.         | Swiss-Prot: Q05519NCBI Gene ID: 9295                                                                                                               |
| Uniprot               | Q05519                                                                                                                                             |
| GeneID                | 9295;                                                                                                                                              |
| SDS-PAGE MW           | 54kd                                                                                                                                               |
| Concentration         | 1.0mg/ml                                                                                                                                           |
| Formulation           | Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage               | Store at -20°C                                                                                                                                     |

## Application Details

Western blotting: 1:500~1:3000

## Images



Western blot analysis of extracts from LOVO cells and RAW264.7 cells, using SFRS11 antibody #35073.

## Background

---

May function in pre-mRNA splicing.

Chaudhary N., Proc. Natl. Acad. Sci. U.S.A. 88:8189-8193(1991).

Gregory S.G., Gregory S.G.,

Page-McCaw P.S., RNA 5:1548-1560(1999).

---

Note: This product is for in vitro research use only