

STK39 Antibody

Catalog No: #35079

Package Size: #35079-1 50ul #35079-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	STK39 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total STK39 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human STK39.
Target Name	STK39
Other Names	DCHT; EC 2.7.11.1; Serine/threonine-protein kinase 39; SPAK; Ste-20 related kinase
Accession No.	Swiss-Prot: Q9UEW8NCBI Gene ID: 27347
Uniprot	Q9UEW8
GeneID	27347;
SDS-PAGE MW	60kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from COLO cells, using STK39 antibody #35079.

Background

May act as a mediator of stress-activated signals.

Johnston A.M., Oncogene 19:4290-4297(2000).

Melnick M.B., Submitted (OCT-1997) to the EMBL/GenBank/DDBJ databases.

Olsen J.V., Cell 127:635-648(2006).

Note: This product is for in vitro research use only