

TBX18 Antibody

Catalog No: #35087

Package Size: #35087-1 50ul #35087-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TBX18 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total TBX18 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human TBX18.
Target Name	TBX18
Other Names	T-box protein 18; T-box transcription factor TBX18; TBX18;
Accession No.	Swiss-Prot: O95935NCBI Gene ID: 9096
Uniprot	O95935
GeneID	9096;
SDS-PAGE MW	68kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from Jurkat cells, using TBX18 antibody #35087.

Background

Probable transcriptional regulator involved in developmental processes.

Mungall A.J., Nature 425:805-811(2003).

The MGC Project Team, Genome Res. 14:2121-2127(2004).

Yi C.-H., Genomics 55:10-20(1999).

Note: This product is for in vitro research use only