

CD6 Antibody

Catalog No: #35088

Package Size: #35088-1 50ul #35088-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CD6 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total CD6 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human CD6.
Target Name	CD6
Other Names	T12; TP120;
Accession No.	Swiss-Prot: P30203NCBI Gene ID: 923
Uniprot	P30203
GeneID	923;
SDS-PAGE MW	73kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from Jurkat cells, using CD6 antibody #35088.

Background

Involved in cell adhesion. Binds to CD166.

Sakaki Y., Nature 440:497-500(2006).

Rudd C.E., J. Biol. Chem. 266:7137-7143(1991)

Aruffo A., J. Exp. Med. 177:219-223(1993)..

Note: This product is for in vitro research use only