

TCEAL3/5/6 Antibody

Catalog No: #35100

Package Size: #35100-1 50ul #35100-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TCEAL3/5/6 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total TCEAL3/5/6 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human TCEAL3/5/6.
Target Name	TCEAL3/5/6
Other Names	Transcription elongation factor A protein-like 3; TCEA-like protein 3; Transcription elongation factor S-II protein-like 3; TCEAL3; MSTP072
Accession No.	Swiss-Prot: Q969E4/Q5H9L2/Q6IPX3NCBI Gene ID: 85012
Uniprot	Q969E4
GeneID	85012;
SDS-PAGE MW	23kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:50~1:100

Images



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue using TCEAL3/5/6 antibody #35100.

Background

May be involved in transcriptional regulation.

Sheng H., Submitted (JUN-1999) to the EMBL/GenBank/DDBJ databases.

Ross M.T., Nat. Genet. 36:40-45(2004).

Ross M.T., Nature 434:325-337(2005).

Note: This product is for in vitro research use only