

TCEAL4 Antibody

Catalog No: #35101

Package Size: #35101-1 50ul #35101-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TCEAL4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total TCEAL4 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human TCEAL4.
Target Name	TCEAL4
Other Names	Transcription elongation factor A protein-like 4; TCEA-like protein 4; Transcription elongation factor S-II protein-like 4; TCEAL4; NPD017
Accession No.	Swiss-Prot: Q96EI5NCBI Gene ID: 79921
Uniprot	Q96EI5
GeneID	79921;
SDS-PAGE MW	25kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:50~1:100

Images



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue using TCEAL4 antibody #35101.

Background

May be involved in transcriptional regulation.

Lu X.W., Zhongguo Yi Xue Ke Xue Yuan Xue Bao 23:27-31(2001).

Xu X., Submitted (MAY-2000) to the EMBL/GenBank/DDBJ databases.

Ota T., Nat. Genet. 36:40-45(2004).

Note: This product is for in vitro research use only