

USP38 Antibody

Catalog No: #35122

Package Size: #35122-1 50ul #35122-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	USP38 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total USP38 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human USP38.
Target Name	USP38
Other Names	Ubiquitin carboxyl-terminal hydrolase 38; Ubiquitin thioesterase 38; Ubiquitin-specific-processing protease 38; Deubiquitinating enzyme 38; HP43.8KD
Accession No.	Swiss-Prot: Q8NB14NCBI Gene ID: 84640
Uniprot	Q8NB14
GeneID	84640;
SDS-PAGE MW	117kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from 293 cells, using USP38 antibody #35122.

Background

Deubiquitinating enzyme exhibiting a preference towards 'Lys-63'-linked Ubiquitin chains.

Ota T., Nat. Genet. 36:40-45(2004).

Bechtel S., BMC Genomics 8:399-399(2007).

Mural R.J., Submitted (SEP-2005) to the EMBL/GenBank/DDBJ databases.

Note: This product is for in vitro research use only