

## USP42 Antibody

Catalog No: #35123

Package Size: #35123-1 50ul #35123-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	USP42 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total USP42 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human USP42.
Target Name	USP42
Other Names	deubiquitinating enzyme 42; EC 3.1.2.15; UB42; ubiquitin carboxyl-terminal hydrolase 42; ubiquitin thiolesterase 42
Accession No.	Swiss-Prot: Q9H9J4NCBI Gene ID: 84132
Uniprot	Q9H9J4
GeneID	84132;
SDS-PAGE MW	150kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:500~1:3000

## Images



Western blot analysis of extracts from HT-29 cells, using USP42 antibody #35123.

## Background

---

Deubiquitinating enzyme which may play an important role during spermatogenesis By similarity.  
Western blot analysis of extracts from HT-29 cells, using USP42 antibody.

---

Note: This product is for in vitro research use only