

UBFD1 Antibody

Catalog No: #35128

Package Size: #35128-1 50ul #35128-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	UBFD1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous levels of total UBFD1 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human UBFD1.
Target Name	UBFD1
Other Names	FLJ38870; FLJ42145; ubiquitin family domain containing 1; UBPH;
Accession No.	Swiss-Prot: O14562NCBI Gene ID: 56061
Uniprot	O14562
GeneID	56061;
SDS-PAGE MW	50kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from K562 cells, using UBFD1 antibody #35128.

Background

May play a role as NF-kappa-B regulator.

Loftus B.J., Genomics 60:295-308(1999).

Ota T., Nat. Genet. 36:40-45(2004).

Note: This product is for in vitro research use only