

ZNF460 Antibody

Catalog No: #35158

Package Size: #35158-1 50ul #35158-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ZNF460 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total ZNF460 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human ZNF460.
Target Name	ZNF460
Other Names	HZF8; zinc finger protein 272; zinc finger protein HZF8; ZN272; ZN460
Accession No.	Swiss-Prot: Q14592NCBI Gene ID: 10794
Uniprot	Q14592
GeneID	10794;
SDS-PAGE MW	64kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from LOVO cells, using ZNF460 antibody #35158.

Background

May be involved in transcriptional regulation.

Dai J., Cytogenet. Genome Res. 103:74-78(2003).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

Abrink M., DNA Cell Biol. 14:125-136(1995).

Note: This product is for in vitro research use only