

SLC30A9 Antibody

Catalog No: #35161

Package Size: #35161-1 50ul #35161-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SLC30A9 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total SLC30A9 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human SLC30A9.
Target Name	SLC30A9
Other Names	Zinc transporter 9ZnT-9; Solute carrier family 30 member 9; Human embryonic lung protein; HUEL SLC30A9;
Accession No.	Swiss-Prot: Q6PML9NCBI Gene ID: 10463
Uniprot	Q6PML9
GeneID	10463;
SDS-PAGE MW	63kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from HepG2 cells, using SLC30A9 antibody #35161.

Background

Plays a role in the p160 coactivator signaling pathway that mediates transcriptional activation by nuclear receptors By similarity. Plays a role in transcriptional activation of Wnt-responsive genes By similarity.

Sim D.L.C., Genomics 59:224-233(1999).

Chow V.T.K., Submitted (APR-2000) to the EMBL/GenBank/DDBJ databases.

Xu Z.Y., Submitted (JUN-2003) to the EMBL/GenBank/DDBJ databases.

Note: This product is for in vitro research use only