

TTF2 Antibody

Catalog No: #35169

Package Size: #35169-1 50ul #35169-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TTF2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IF
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total TTF2 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human TTF2.
Target Name	TTF2
Other Names	Interferon-related developmental regulator 1; Nerve growth factor-inducible protein PC4;
Accession No.	Swiss-Prot: O00358NCBI Gene ID: 2304
Uniprot	O00358
GeneID	2304;
SDS-PAGE MW	50kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Immunofluorescence: 1:100~1:500

Images



Immunofluorescence analysis of A549 cells, using TTF2 antibody #35169.

Background

Probable transcription factor. Could be involved in thyroid gland organogenesis.

Buane P., Genomics 51:233-242(1998).

Hillier L.W., Nature 424:157-164(2003).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

Note: This product is for in vitro research use only