

ZNF498 Antibody

Catalog No: #35172

Package Size: #35172-1 50ul #35172-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ZNF498 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total ZNF498 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human ZNF498.
Target Name	ZNF498
Other Names	Zinc finger protein 498; Zinc finger and SCAN domain-containing protein 25; ZNF498; ZSCAN25;
Accession No.	Swiss-Prot: Q6NSZ9NCBI Gene ID: 221785
Uniprot	Q6NSZ9
GeneID	221785;
SDS-PAGE MW	30kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from COS-7 cells, using ZNF498 antibody #35172.

Background

May be involved in transcriptional regulation By similarity.

Scherer S.W., Science 300:767-772(2003).

Mural R.J., Submitted (JUL-2005) to the EMBL/GenBank/DDBJ databases.

The MGC Project Team; Genome Res. 14:2121-2127(2004).

Note: This product is for in vitro research use only