

RALY Antibody

Catalog No: #35179

Package Size: #35179-1 50ul #35179-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RALY Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total RALY protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human RALY.
Target Name	RALY
Other Names	autoantigen p542; hnRNP associated with lethal yellow homolog; P542; RNA-binding protein Raly;
Accession No.	Swiss-Prot: Q9UKM9NCBI Gene ID: 22913
Uniprot	Q9UKM9
GeneID	22913;
SDS-PAGE MW	35kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from Jurkat cells and COLO205 cells, using RALY antibody #35179.

Background

Probable-RNA binding protein. Could be a heterogeneous nuclear ribonucleoprotein (hnRNP). May be involved in pre-mRNA splicing By similarity.

Khrebtukova I., Biochim. Biophys. Acta 1447:107-112(1999).

Rhodes G.H., J. Autoimmun. 10:447-454(1997).

Deloukas P., Nature 414:865-871(2001).

Note: This product is for in vitro research use only