

A20A1 Antibody

Catalog No: #35180

Package Size: #35180-1 50ul #35180-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	A20A1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total A20A1 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human A20A1.
Target Name	A20A1
Other Names	A20A1; ANKRD20A; ANKRD20A1; ankyrin repeat domain 20 family; member A1
Accession No.	Swiss-Prot: Q5TYW2NCBI Gene ID: 84210/441430/441425/728747
Uniprot	Q5TYW2
GeneID	84210;
SDS-PAGE MW	95kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from 3T3 cells, using A20A1 antibody #35180.

Background

Wiemann S., Genome Res. 11:422-435(2001) [PubMed: 11230166].

Humphray S.J., Nature 429:369-374(2004) [PubMed: 15164053].

Mayya V., Sci. Signal. 2:RA46-RA46(2009) [PubMed: 19690332].

Note: This product is for in vitro research use only