

ZNF691 Antibody

Catalog No: #35181

Package Size: #35181-1 50ul #35181-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ZNF691 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total ZNF691 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from Internal of human ZNF691.
Target Name	ZNF691
Other Names	Zinc finger protein 691; ZNF691;
Accession No.	Swiss-Prot: Q5VV52NCBI Gene ID: 51058
Uniprot	Q5VV52
GeneID	51058;
SDS-PAGE MW	36kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from HeLa cells, using ZNF691 antibody #35181.

Background

May be involved in transcriptional regulation.

Mei G., Submitted (FEB-1999) to the EMBL/GenBank/DDBJ databases.

Ota T., Nat. Genet. 36:40-45(2004).

Ebert L., Submitted (JUN-2004) to the EMBL/GenBank/DDBJ databases.

Note: This product is for in vitro research use only