

ZNF446 Antibody

Catalog No: #35203

Package Size: #35203-1 50ul #35203-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ZNF446 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total ZNF446 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human ZNF446.
Target Name	ZNF446
Other Names	Zinc finger protein 446; Zinc finger protein with KRAB and SCAN domains 20; ZNF446; ZKSCAN20;
Accession No.	Swiss-Prot: Q9NWS9NCBI Gene ID: 55663
Uniprot	Q9NWS9
GeneID	55663;
SDS-PAGE MW	49kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from HuvEc cells, using ZNF446 antibody #35203.

Background

May be involved in transcriptional regulation.

Yu L., Submitted (APR-2003).

Sugano S., Nat. Genet. 36:40-45(2004).

Han D.K., Sci. Signal. 2:RA46-RA46(2009)

Note: This product is for in vitro research use only