

FGFRL1 Antibody

Catalog No: #35217

Package Size: #35217-1 50ul #35217-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	FGFRL1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total FGFRL1 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from Internal of human FGFRL1.
Target Name	FGFRL1
Other Names	FGF homologous factor receptor; FGF receptor-like protein 1; FGFR-5; FGFR-like protein; FGFR5
Accession No.	Swiss-Prot: Q8N441NCBI Gene ID: 53834
Uniprot	Q8N441
GeneID	53834;
SDS-PAGE MW	55kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from NIH-3T3 cells, using FGFRL1 antibody #35217.

Background

Has a negative effect on cell proliferation By similarity.

Wiedemann M., Genomics 69:275-279(2000).

Kim I., Biochim. Biophys. Acta 1518:152-156(2001).

Aggarwal S., Submitted (OCT-2000) to the EMBL/GenBank/DDBJ databases.

Note: This product is for in vitro research use only