

TMEM101 Antibody

Catalog No: #35251

Package Size: #35251-1 50ul #35251-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TMEM101 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total TMEM101 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human TMEM101.
Target Name	TMEM101
Other Names	Transmembrane protein 101; Putative NF-kappa-B-activating protein 130; TMEM101;
Accession No.	Swiss-Prot: Q96IK0NCBI Gene ID: 84336
Uniprot	Q96IK0
GeneID	84336;
SDS-PAGE MW	29kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from HuvEc cells, using TMEM101 antibody #35251.

Background

May activate NF-kappa-B signaling pathways.

Matsuda A., Oncogene 22:3307-3318(2003).

Ota T., Nat. Genet. 36:40-45(2004)

Note: This product is for in vitro research use only