

LAMP1 Antibody

Catalog No: #35308



Package Size: #35308-1 50ul #35308-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	LAMP1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total LAMP1 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human LAMP1.
Target Name	LAMP1
Other Names	Lysosome-associated membrane glycoprotein 1; Lysosome-associated membrane protein 1; LAMP-1; CD107 antigen-like family member A; CD107a
Accession No.	Swiss-Prot: P11279NCBI Gene ID: 3916
Uniprot	P11279
GeneID	3916;
SDS-PAGE MW	125kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from K562, COLO205, HeLa, Jurkat and HepG2 cells, using LAMP1 antibody #35308.

Background

Presents carbohydrate ligands to selectins. Also implicated in tumor cell metastasis.

ss M.T., Nature 428:522-528(2004).

Venter J.C., Submitted (JUL-2005).

Hasilik A., Traffic 8:1676-1686(2007).

Note: This product is for in vitro research use only