

FAM3D Antibody

Catalog No: #35324

Package Size: #35324-1 50ul #35324-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	FAM3D Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total FAM3Dprotein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from C-terminal of human FAM3D.
Target Name	FAM3D
Other Names	cytokine-like protein EF-7; EF7; FAM3D; family with sequence similarity 3; member D
Accession No.	Swiss-Prot: Q96BQ1NCBI Gene ID: 131177
Uniprot	Q96BQ1
GeneID	131177;
SDS-PAGE MW	25kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:3000

Images



Western blot analysis of extracts from K562 cells, using FAM3D antibody #35324.

Background

Zhu Y., Genomics 80:144-150(2002).

Clark H.F., Genome Res. 13:2265-2270(2003).

Ramachandran P., J. Proteome Res. 5:1493-1503(2006).

Note: This product is for in vitro research use only