

BCAT2 antibody

Catalog No: #22010



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


Description

Product Name	BCAT2 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC
Species Reactivity	Hu Ms
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 1 and 302 of BCAT2
Target Name	BCAT2
Accession No.	Swiss-Prot:O15382Gene ID:587
Uniprot	O15382
GeneID	587;
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 44kd
Western blotting: 1:500-1:3000
Immunohistochemistry: 1:100-1:250
Immunofluorescence: 1:100-1:200

Images



Sample (30 ug of whole cell lysate)  
A: Raji  
10% SDS PAGE  
Primary antibody diluted at 1: 1000



A:Immunohistochemical analysis of paraffin-embedded Ca922 xenograft, using BCAT2 antibody at 1: 100 dilution.

B:Immunohistochemical analysis of paraffin-embedded Ca922 xenograft, using BCAT2 antibody at 1: 500 dilution.

---

Note: This product is for in vitro research use only