

NCAM Antibody

Catalog No: #35377

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NCAM Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Rt Ms
Specificity	The antibody detects endogenous levels of total NCAM protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fragment corresponding to a region within amino acids 1 and 171 of NCAM.
Target Name	NCAM
Other Names	neural cell adhesion molecule 1 antibody; Cd56 antibody; MGC124601 antibody; N-CAM antibody; N-CAM-1 antibody; NCAM-1 antibody; NCAM-C antibody; NCAMC antibody; Ncam antibody; Ncam1 antibody
Accession No.	Swiss-Prot#:P13596;NCBI Gene#:24586
Uniprot	P13596
GeneID	24586;
SDS-PAGE MW	95kd
Concentration	1mg/ml
Formulation	Rabbit IgG in 1XPBS, 40% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C


Application Details

Western blotting: 1:500-1:3000

Images



Sample (50 ug of whole cell lysate)
A: Rat brain
5% SDS PAGE
#35377 diluted at 1:2000



Sample (50 ug of whole cell lysate)

A: mouse brain

5% SDS PAGE

#35377 diluted at 1:2000

Background

cell adhesion molecule; involved in neuronal cell adhesion [RGD]

Note: This product is for in vitro research use only