

FBXW11 Antibody

Catalog No: #35405



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	FBXW11 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total FBXW11 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human CENPC (NP_001803.2).
Target Name	FBXW11
Other Names	CENPC;CENP-C;CENPC1;MIF2;hcp-4
Accession No.	Uniprot:Q03188GeneID:1060
Uniprot	Q03188
GeneID	1060
SDS-PAGE MW	130kDa
Concentration	0.19mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using CENPC antibody.

Background

Centromere protein C 1 is a centromere autoantigen and a component of the inner kinetochore plate. The protein is required for maintaining proper

kinetochore size and a timely transition to anaphase. A putative pseudogene exists on chromosome 12.

Note: This product is for in vitro research use only