

FLIP Antibody

Catalog No: #35423



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	FLIP Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total FLIP protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide corresponding to a region within amino acids 418 and 480 of FLIP.
Target Name	FLIP
Other Names	CASH antibody; CASP8AP1 antibody; CLARP antibody; Casper antibody; FLAME antibody; FLAME-1 antibody; FLAME1 antibody; FLIP antibody; I-FLICE antibody; MRIT antibody; c-FLIP antibody; c-FLIPL antibody; c-FLIPR antibody; c-FLIPS antibody; CFLAR antibody; ca
Accession No.	Swiss-Prot#:O15519;NCBI Gene#:8837
Uniprot	O15519
GeneID	8837;
SDS-PAGE MW	55kd
Concentration	1mg/ml
Formulation	Rabbit IgG in 1XPBS, 40% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C

Application Details

Western blotting: 1:500-1:3000

Images



Sample (30 ug of whole cell lysate)
A: A431
7.5% SDS PAGE
#35423 diluted at 1:1000

Background

Apoptosis regulator protein which may function as a crucial link between cell survival and cell death pathways in mammalian cells. Acts as an inhibitor of TNFRSF6 mediated apoptosis. A proteolytic fragment (p43) is likely retained in the death-inducing signaling complex (DISC) thereby blocking further recruitment and processing of caspase-8 at the complex. Full length and shorter isoforms have been shown either to induce apoptosis or to reduce TNFRSF-triggered apoptosis. Lacks enzymatic (caspase) activity.

Note: This product is for in vitro research use only