

FRK Antibody

Catalog No: #35488



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	FRK Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total FRK protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fragment contain a sequence corresponding to a region within amino acids 226 and 430 of FRK.
Target Name	FRK
Other Names	GTK antibody; PTK5 antibody; RAK antibody; FRK antibody; tyrosine-protein kinase FRK antibody; nuclear tyrosine protein kinase RAK antibody; protein-tyrosine kinase 5 antibody; PTK5 protein tyrosine kinase 5 antibody; fyn-related kinase antibody
Accession No.	Swiss-Prot#:P42685;NCBI Gene#:2444
Uniprot	P42685
GeneID	2444;
SDS-PAGE MW	58kd
Concentration	0.98mg/ml
Formulation	Rabbit IgG in 0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C

Application Details

Western blotting: 1:500-1:3000

Images



Sample (30 ug of whole cell lysate)
A: Hep G2
7.5% SDS PAGE
#35488 diluted at 1:1000

Background

The protein encoded by this gene belongs to the TYR family of protein kinases. This tyrosine kinase is a nuclear protein and may function during G1 and S phase of the cell cycle and suppress growth. [provided by RefSeq]

Note: This product is for in vitro research use only