

UPK3B Antibody

Catalog No: #35546



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description


Product Name	UPK3B Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	WB IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total UPK3B protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	UPK3B
Other Names	FLJ32198; MGC10902; p35; UPIIIb
Accession No.	Swiss-Prot#: Q9BT76NCBI Gene ID: 80761Gene Accssion: BC004304
Uniprot	Q9BT76
GeneID	105375355;
SDS-PAGE MW	34kd
Concentration	1.1mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN3, 40% Glycerol.
Storage	Store at -20°C

Application Details


Western blotting: 1:1000-1:4000

Immunohistochemistry: 1:25-1:50

Images



Gel: 10%SDS-PAGE
Lysate: 40ug 293T cell
Primary antibody: 1/1000 dilution
Secondary antibody dilution: 1/8000
Exposure time: 1 minute



Immunohistochemical analysis of paraffin-embedded Human lung cancer tissue using #35546 at dilution 1/10.

Background

UPK3B is a minor component of the apical plaques of mammalian urothelium that binds and dimerizes with uroplakin-1b (UPK1B, MIM 602380), one of the major conserved urothelium membrane proteins. The other major conserved integral membrane proteins of urothelial plaques are UPK1A (MIM 611557), UPK2 (MIM 611558), and UPK3A (MIM 611559)

Note: This product is for in vitro research use only