

IL17RD Antibody

Catalog No: #35602

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	IL17RD Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total IL17RD protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from internal residues of human Interleukin-17 receptor D
Target Name	IL17RD
Other Names	SEF; IL-17RD; IL17RLM
Accession No.	Swiss-Prot#: Q8NFM7NCBI Gene ID: 54756Gene Accssion: BC111702
Uniprot	Q8NFM7
GeneID	54756;
Concentration	0.2mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:5-1:20

Images



Immunohistochemical analysis of paraffin-embedded Human lung cancer tissue using #35602 at dilution 1/5.

Background

This gene encodes a membrane protein belonging to the interleukin-17 receptor (IL-17R) protein family. The encoded protein is a component of the interleukin-17 receptor signaling complex, and the interaction between this protein and IL-17R does not require the interleukin (PMID: 19079364). The gene product also affects fibroblast growth factor signaling, inhibiting or stimulating growth through MAPK/ERK signaling.

Note: This product is for in vitro research use only