

ABCB5 Antibody

Catalog No: #35612

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ABCB5 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total ABCB5 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from internal residues of human ATP-binding cassette sub-family B member 5
Target Name	ABCB5
Other Names	ABCB5beta; EST422562; ABCB5alpha
Accession No.	Swiss-Prot#: Q2M3G0NCBI Gene ID: 340273Gene Accssion: BC104894
Uniprot	Q2M3G0
GeneID	340273;
Concentration	0.6mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN3, 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1:10-1:50

Images



Immunohistochemical analysis of paraffin-embedded Human cervical cancer tissue using #35612 at dilution 1/10.



Immunohistochemical analysis of paraffin-embedded Human renal cancer tissue using #35612 at dilution 1/10.

Background

ABCB5 belongs to the ATP-binding cassette (ABC) transporter superfamily of integral membrane proteins. These proteins participate in ATP-dependent transmembrane transport of structurally diverse molecules ranging from small ions, sugars, and peptides to more complex organic molecules.

Note: This product is for in vitro research use only