

## ANXA3 Antibody

Catalog No: #35619

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


Product Name	ANXA3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	WB IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total ANXA3 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Fusion protein corresponding to residues near the N terminal of human annexin A3
Target Name	ANXA3
Other Names	ANX3
Accession No.	Swiss-Prot#: P12429NCBI Gene ID: 306Gene Accssion: BC000871
Uniprot	P12429
GeneID	306;
SDS-PAGE MW	36kd
Concentration	1.3mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.
Storage	Store at -20°C

## Application Details


Western blotting: 1:500-1:2000

Immunohistochemistry: 1:50-1:200

## Images



Gel: 10%SDS-PAGE  
Lysate: 40ug Lovo cell  
Primary antibody: 1/650 dilution  
Secondary antibody dilution: 1/8000  
Exposure time: 1 second



Immunohistochemical analysis of paraffin-embedded Human thyroid cancer tissue using #35619 at dilution 1/35.

## Background

This gene encodes a member of the annexin family. Members of this calcium-dependent phospholipid-binding protein family play a role in the regulation of cellular growth and in signal transduction pathways. This protein functions in the inhibition of phospholipase A2 and cleavage of inositol 1,2-cyclic phosphate to form inositol 1-phosphate. This protein may also play a role in anti-coagulation.

Note: This product is for in vitro research use only